Size:420x280mm

SARS-CoV-2 & Influenza A+B Antigen Combo Rapid Test Cassette Package Insert

Specimens: Nasopharyngeal Swab/
Oropharyngeal Swab
Effective Date: 2022.12
For professional *in vitro* diagnostic use only.

Cat: COIF-522

Version: B

Code: 4.15.03.90003-2

INTENDED USE

The SARS-CoV-2 & Influenza A+B Antigen Combo Rapid Test Cassette is a lateral flow immunoassay intended for the qualitative detection of SARS-CoV-2, influenza A and influenza B viral nucleoprotein antigens from individuals suspected of respiratory viral infection consistent with SARS-CoV-2 by their healthcare provider. Symptoms of Respiratory viral infection due to SARS-CoV-2 and influenza can be similar. The SARS-CoV-2 & Influenza A+B Antigen Combo Rapid Test Cassette is intended for the detection and differentiation of SARS-CoV-2, influenza A and influenza B viral nucleoprotein antigens. Antigens are generally detectable in nasopharyngeal/oropharyngeal specimens during the acute phase of infection. Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. Negative results do not rule out SARS-CoV-2 or influenza infection and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results must be combined with clinical observations, patient history and epidemiological information, and confirmed with a molecular assay, if necessary, for patient management. The SARS-CoV-2 & Influenza A+B Antigen Combo Rapid Test Cassette is intended for use by trained clinical laboratory personnel specifically instructed and trained in *in vitro* diagnostic procedures.

SUMMARY

The novel corona viruses (SARS-CoV-2) belong to the β genus. SARS-CoV-2 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel corona virus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, malignant diarrhea is found in a few cases. Influenza (Flu) is a contagious respiratory illness caused by influenza viruses. It can cause mild to severe illness. Serious outcomes of flu infection can result in hospitalization, even death, especially among children, and people with certain health conditions are at high risk of serious flu complications. There are two main types of influenza (flu) virus: Types A and B. The influenza A and viruses that routinely spread in people (human influenza viruses) are responsible for seasonal flu epidemics each year.

DETAILS

The SARS-CoV-2 & Influenza A+B Antigen Combo Rapid Test Cassette is a qualitative, lateral flow immunoassay for the detection of the N protein of SARS-CoV-2, Influenza A and Influenza B nucleoproteins in Nasopharyngeal swab/oropharyngeal swab. In this test, antibody specific to the N protein of SARS-CoV-2, Influenza A and Influenza B nucleoproteins is separately coated on the test line regions of the test cassette. During testing, the extracted specimen reacts with the antibody to N protein of SARS-CoV-2, Influenza A and/or Influenza B that are coated onto particles. The mixture migrates up the membrane to react with the antibody to N protein of SARS-CoV-2, Influenza A and/or Influenza B on the membrane and generate one colored line in the test regions. The presence of this colored line of the test regions indicates a positive result. To serve as a procedural control, a colored line will always appear in the control region if the test has performed properly.

REAGENTS

The test cassette contains anti-SARS-CoV-2 Nucleocapsid, anti-Influenza A and B protein particles, and anti-SARS-CoV-2 Nucleocapsid, anti-Influenza A and B protein coated on the membrane.

PRECAUTIONS

Please read all the information in this package insert before performing the test.

- For professional *in vitro* diagnostic use only. Do not use after the expiration date.
- The test should remain in the sealed pouch until ready to use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infection agent.
- The used test should be discarded according to local regulations.
- Avoid using bloody samples.
- Wear gloves when handling the samples, avoid touching the reagent membrane and sample well.

STORAGE AND STABILITY

1. Store as packaged in the sealed pouch at the temperature (2-30°C). The kit is stable within the expiration date printed on the labeling. **DO NOT FREEZE**.

2. Once open the pouch, the test should be used within 15 min. Prolonged exposure to hot and humid environment will cause product deterioration. The lot number and the expiration date were printed on the labeling.

SPECIMEN COLLECTION, TRANSPORT AND STORAGE

Specimen Collection:

Nasal secretion collection:
1. Insert a sterile swab into the nostril of the patient, reaching the surface of the posterior nasopharynx.
2. Swab over the surface of the posteriornasopharynx.
3. Withdraw the sterile swab from the nasal cavity.

Throat secretion collection:

1. Insert a sterile swab into the throat completely from the mouth, centering on the throat wall and the reddened area of the palatal tonsils, wipe the bilateral pharyngeal tonsils and posterior pharyngeal wall with moderate force, avoid touching the tongue and take out the swab.

Specimen transport and storage:

- Specimens should be tested as soon as possible after collection.
- If swabs are not processed immediately, it is highly recommended the swab sample is placed into a dry, sterile, and tightly sealed plastic tube for storage. It can be stored at 2-8°C for 8 hours and can be stored for a long time at -70°C. Do not store specimens in viral transport media.
- Samples collected from swabs that are too viscous or agglomerated are not recommended for testing of this product. If the swabs are contaminated with a large amount of blood, they are not recommended for testing. It is not recommended to use the samples that are processed with sample extraction solution not provided in this kit for testing of this product.

SPECIMEN PREPARATION

Only the extraction buffer and tubes provided in the kit is to be used for swab specimen preparation.

- Take out the tube, peel off the aluminum foil on the extraction buffer tube carefully avoid liquid spills and place the tube in the Tube Holder.
- Place the swab specimen in the Extraction Tube. Rotate the swab for approximately 10 seconds while pressing the head against the inside of the tube to release the antigen in the swab.
- Remove the swab while squeezing the swab head against the inside of the Extraction Tube as you remove it to expel as much liquid as possible form the swab. Discard the swab in accordance with your biohazard waste disposal protocol.
- Fit the dropper tip on top of the extraction tube. Place the test cassette on a clean and level surface.

INTERPRETATION OF RESULTS
(Please refer to the illustration above)

- POSITIVE SARS-CoV-2+Two distinct colored lines appear in the left window.**
One colored line should be in the control region (C) and another colored line should be in the Test region (T). Positive result in the Test region indicates detection of SARS-CoV-2 antigens in the sample.
- POSITIVE Influenza A+Two distinct colored lines appear in the right window.**
One colored line should be in the control region (C) and another colored line should be in the Influenza A region (A). A positive result in the Influenza A region indicates that Influenza A antigen was detected in the sample.
- POSITIVE Influenza B+Two distinct colored lines appear in the right window.**
One colored line should be in the control region (C) and another colored line should be in the Influenza B region (B). A positive result in the Influenza B region indicates that Influenza B antigen was detected in the sample.

- SPECIMEN PREPARATION**
The test cassette contains anti-SARS-CoV-2, Influenza A and/or Influenza B that are coated onto particles. The mixture migrates up the membrane to react with the antibody to N protein of SARS-CoV-2, Influenza A and/or Influenza B on the membrane and generate one colored line in the test regions. The presence of this colored line of the test regions indicates a positive result. To serve as a procedural control, a colored line will always appear in the control region if the test has performed properly.

REAGENTS
The test cassette contains anti-SARS-CoV-2 Nucleocapsid, anti-Influenza A and B protein particles, and anti-SARS-CoV-2 Nucleocapsid, anti-Influenza A and B protein coated on the membrane.

PRECAUTIONS
The test readout and the test should be considered with other clinical findings from other laboratory tests and evaluations.

LIMITATIONS OF THE TEST
1. The test Procedure and the Interpretation of test Result must be followed closely when testing for the presence of SARS-CoV-2/Influenza A/Influenza B antigens in the human nasopharyngeal/oropharyngeal specimens from suspected individuals. For optimal test performance, proper sample collection is critical. Failure to follow the procedure may give inaccurate results.

2. The performance of the SARS-CoV-2 & Influenza A+B Antigen Combo Rapid Test Cassette was evaluated using the procedures provided in this product insert only. Modifications to these procedures may alter the performance of the test. Viral Transport Media (VTM) specimen and extracted specimens for PCR tests cannot be used for the test.

3. The SARS-CoV-2 & Influenza A+B Antigen Combo Rapid Test Cassette is for *in vitro* diagnostic use only. This test should be used for detection of SARS-CoV-2/Influenza A/B Antigens in human nasopharyngeal/oropharyngeal specimens as an aid in the diagnosis of patients with suspected SARS-CoV-2, Influenza A or Influenza B infection in conjunction with clinical presentation and the results of other laboratory tests. Neither the quantitative value nor the rate of increase in the concentration of SARS-CoV-2/Influenza A/Influenza B antigens can be determined by this qualitative test.

4. The SARS-CoV-2 & Influenza A+B Antigen Combo Rapid Test Cassette will only indicate the presence of SARS-CoV-2/Influenza A/Influenza B Antigens in the specimen and should not be used as the sole criteria for the diagnosis of SARS-CoV-2/Influenza A/Influenza B infection.

5. Test results obtained with the test should be considered with other clinical findings from other laboratory tests and evaluations.

6. If the test result is negative, and clinical symptoms persist, it is recommended to re-sample the patient a few days later and test again or test with a molecular diagnostic device to rule out infection in these individuals.

7. The test will show negative results under the following conditions: The concentration of the novel coronavirus antigens, Influenza A or Influenza B virus in the sample is lower than the minimum detection limit of the test.

8. Negative results do not rule out SARS-CoV-2 infection, particularly in those who have been in contact with the virus. Follow-up testing with a molecular diagnostic should be considered to rule out infection in these individuals.

9. A negative result for Influenza A or Influenza B obtained from this kit should be confirmed by RT-PCR culture.

10. Excess blood or mucin on the swab specimen may interfere with test performance and may yield a false positive result.

11. The accuracy of the test depends on the quality of the swab sample. False negatives may result from improper sample collection or storage.

12. Positive results of SARS-CoV-2 may be due to infection with non-SARS-CoV-2 coronavirus strains or other interference factors. A positive result for influenza A and/or B does not preclude an underlying co-infection with another pathogen, therefore the possibility of an underlying bacterial infection should be considered.

EXPECTED VALUES

The SARS-CoV-2 & Influenza A+B Antigen Combo Rapid Test Cassette has been compared with leading commercial RT-PCR tests. The correlation between these two systems is no less than 97%.

PERFORMANCE CHARACTERISTICS

Sensitivity, Specificity and Accuracy

The SARS-CoV-2 & Influenza A+B Antigen Combo Rapid Test Cassette has been evaluated with specimens obtained from the patients' PCR is used as the reference method for the SARS-CoV-2 & Influenza A+B Antigen Combo Rapid Test Cassette. Specimens were considered positive if RT-PCR indicated a positive result. Specimens were considered negative if RT-PCR indicated a negative result.

SPECIMEN COLLECTION, TRANSPORT AND STORAGE

Description Cross reaction

Human adenovirus 3 N/A

Human adenovirus 7 N/A

Human coronavirus OC43 N/A

Human metapneumovirus N/A

Human rhinovirus N/A

Human respiratory syncytial virus N/A

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