

Glucosity PC Software



Overview

You can view the overall statistics in the selected time period.

- 7, 14, 30, 90-day averaging
- Low / high glucose marker

Meal and bedtime

You can view statistics of certain periods of a day over several days.



Trend

You can view the tendency of blood glucose over several days.



BG Table

You can check the blood glucose test records of individual days. Values out of the target range is color-marked.

low ● normal ● high ●

- Blood glucose numerical record
- Test value color-coding
- Daily and weekly review



Standard day

Blood glucose values of the selected periods and mapped on a 24-hour scale, you can check the dispersion and discrepancies in an amplitude view.